


 **Movora**
EXCLUSIVE



DemeTECH[®]

"A Trusted Name in Medicine"

Movora is the proud distributor of wound closure products from DemeTECH[®], a leading North American company that produces suture.



Proudly Made in the USA







Table of Contents

Introduction	4
DemeTRACK	5
Needle Options	6
Suture Length and Material	7
Sutures	8-25
Tensile Strength Comparison	9
DemeSORB™	10
DemeDIOX™	12
DemeGUT™	14
DemeCAPRONE™	16
DemeCRYL™	18
DemeLON™	20
DemeLENE™	22
DemeSILK™	24
Antibacterial Sutures	26
Non-Suture Offerings	27

Introduction

DemeTECH Corporation is a world-renowned leader in surgical sutures and blades. DemeTECH strives to enhance doctor-patient relationships through hard work, dedication, and a commitment to the pursuance of cutting-edge technology and innovation. Their corporate headquarters are located in Miami, Florida, with office facilities throughout North America, South America, Asia, the Middle East, and Europe.

DemeTECH is an exceptional corporation, in that it uses next-generation technology to meet the growing demands of today's society. Our quality has been proven time and time again: high demand and worldwide use, distribution in over 100 countries worldwide, and certification of the three major medical authorities: FDA, CE, and ISO.



"So far so good! We have been using the VOI lines of suture and have had zero problems. Handling and functionality seem comparable to more expensive brands. Thanks VOI!"

*Dr. Levien, Australia
Using PDO, PGCL, PGA, Nylon*




DemeTRACK Benefits



DemeTRACK reduces the memory effect observed on monofilament sutures.



DemeTECH's DemeTRACK is available for:

- Monofilament Sutures
- Multifilament Sutures
- Sutures less than 150cm
- Single Thread Sutures
- Double Needle Sutures




The photo illustration on the left compares the reduced memory on a standard suture packaged with DemeTECH'S DemeTRACK and an equivalent suture in conventional packaging.


How to Use DemeTRACK



Needle Length




Needle Type and Profile




Suture Lengths

Purpose	Suture Size
Delicate surgeries	11-0 thru 8-0
Small vessels, facial surgeries or graft sewing	7-0, 6-0
Vessel repair or skin closure	5-0, 4-0
Skin closure, closure of muscle layers	3-0, 2-0
Fascia layer in abdominal surgery, knee and hip surgery	0, 1
Repair of tendons or other high tension structure	2 thru 5


Suture Materials



Sutures





Tensile Strength Comparison In Absorbable Sutures



* Absorption rates are over a period of 7, 14, 21, 28, 42, 56-70, 90-120 and 180 days where applicable between Demetech and the Leading Brand suture types for PGA (Polyglycolic Acid), PDO (Polydioxanone) and PGCL (Poliglecaprone 25.)

DemeSORB™

Polyglycolic Acid Suture




DemeTech's PGA Sutures

DemeSORB Polyglycolic Acid Suture is a synthetic, absorbable, sterile, surgical suture composed of 100% glycolide. The chemical formula of DemeTech PGA is $(C_2H_2O_2)_n$. Clinical trials have shown that after two weeks, approximately 84% of DemeTech's PGA initial strength remained. At four weeks, approximately 23% of DemeTech's PGR linear tensile strength remained. DemeSORB™ Polyglycolic Acid Suture has been found to be nonantigenic and elicit only mild tissue reactivity during absorption process.

Characteristics of DemeTech's PGA Sutures

- Polyglycolic Acid, better known as PGA, is composed of polymers of glycolide acids.
- DemeTech's PGA is coated to allow for easier passage through tissues with only minimal drag.
- Coated PGA Sutures facilitate:
 - Easier Tissue Passage
 - Precise Knot Placement
 - Smooth Tie Down
 - A decreased Tendency to Irritate Tissues
 - Greater Tensile Strength



Tensile Strength

Multifilament



Violet

Colorless



Material

Complete Absorption

60 ► 90 Days



Absorption Profile

Polyglycolic Acid Suture

ITEM #	DESCRIPTION
PGA1082036B0P	DemeSORB, Polyglycolic Acid, 90cm Long, Violet, 2/0, 36mm Needle, 1/2 Circle, RB, DCT-1
PGA1172026B0P	DemeSORB, Polyglycolic Acid, 70cm Long, Colorless, 2/0, 26mm Needle, 1/2 Circle, RB, DSH
PGA1173019F13M	DemeSORB, Polyglycolic Acid, 70cm L, Colorless, 3/0, 19mm Needle, 3/8 Circle Precision Pt RC, DPS-2
PGA1182026B0P	DemeSORB, Polyglycolic Acid, 70cm Long, Violet, 2/0, 26mm Needle, 1/2 Circle, RB, DSH
PGA1183026H0P	DemeSORB, Polyglycolic Acid, 70cm Long, Violet, 3/0, 26mm Needle, DUR-6, 5/8 Circle, RB
PGA1670265B0P	DemeSORB, Polyglycolic Acid, 135cm Long, Colorless, 2, 65mm Needle, 1/2 Circle, RB, DTP-1
PGA176013F13M	DemeSORB, Polyglycolic Acid, 45cm Long, Colorless, 6/0, 13mm Needle, 3/8 Circle, Precision Point RC
PGA183019F4P	DemeSORB, Polyglycolic Acid, 45cm Long, Violet, 3/0, 19mm Needle, 3/8 Circle, RC, DFS-2
PGA184013F13M	DemeSORB Polyglycolic Acid, 45cm Long, Violet, 4/0, 13mm Needle, 3/8 Circle Precision Point RC, DP-3
PGA184013F4P	DemeSORB, Polyglycolic Acid, 45cm Long, Violet, 4/0, 13mm Needle, 3/8 Circle, RC, DC-3
PGA184019F4P	DemeSORB, Polyglycolic Acid, 45cm Long, Violet, 4/0, 19mm Needle, 3/8 Circle, RC, DFS-2
PGA185013F4P	DemeSORB, Polyglycolic Acid, 45cm Long, Violet, 5/0, 13mm Needle, 3/8 Circle, RC, DC-3
PGA185019F4P	DemeSORB, Polyglycolic Acid, 45cm Long, Violet, 5/0, 19mm Needle, 3/8 Circle, RC, DFS-2
PGA186013F13M	DemeSORB Polyglycolic Acid, 45cm Long, Violet, 6/0, 13mm Needle, 3/8 Circle Precision Point RC, DP-3
PGA270126B4P	DemeSORB, Polyglycolic Acid, 75cm Long, Colorless, 1, 26mm, DCP-2, 1/2 Circle, RC
PGA276017B0P	DemeSORB, Polyglycolic Acid, 75cm Long, Colorless, 6/0, 17mm Needle, DRB-1, 1/2 Circle, RB
PGA280026B1P	DemeSORB, Polyglycolic Acid, 75cm Long, Violet, 0, 26mm Needle, DCT-2, 1/2 Circle, RB Heavy
PGA280030B4P	DemeSORB, Polyglycolic Acid, 75cm Long, Violet, 0, 30 MM Needle, 1/2 Circle, RC DCPA
PGA280030F13M	DemeSORB, Polyglycolic Acid, 75cm Long, Violet, 0, 30mm Needle, 3/8 Circle Precision Point RC, DPSL
PGA280036B4P	DemeSORB, Polyglycolic Acid, 75cm Long, Violet, 0, 36mm Needle, 1/2 Circle, RC, DCP-1
PGA280037B0P	DemeSORB, Polyglycolic Acid, 75cm Long, Violet, 0, 37mm Needle, 1/2 Circle, RB, DCT-1
PGA280136B0P	DemeSORB, Polyglycolic Acid, 75cm Long, Violet, 1, 36mm Needle, DCT-1, 1/2 Circle, RB
PGA280137B0P	DemeSORB, Polyglycolic Acid, 75cm Long, Violet, 1, 37mm Needle, 1/2 Circle, RB, DCT-1
PGA282022B4P	DemeSORB, Polyglycolic Acid, 75cm Long, Violet, 2/0, 22mm Needle, DX-1, 1/2 Circle, RC
PGA282024F4P	DemeSORB, Polyglycolic Acid, 75cm Long, Violet, 2/0, 24mm Needle, 3/8 Circle, RC, DFS-1
PGA282026B0P	DemeSORB, Polyglycolic Acid, 75cm Long, Violet, 2/0, 26mm Needle DSH, 1/2 Circle, RB
PGA282026B1P	DemeSORB, Polyglycolic Acid, 75cm Long, Violet, 2/0, 26mm Needle, DCT-2, 1/2 Circle, RB Heavy
PGA282030F13M	DemeSORB Polyglycolic Acid, 75cm Long, Violet, 2/0, 30mm Needle, 3/8 Circle Precision Point RC, DPSL
PGA282036B4P	DemeSORB, Polyglycolic Acid, 75cm Long, Violet, 2/0, 36mm Needle, DCP-1, 1/2 Circle, RC
PGA282037B0P	DemeSORB, Polyglycolic Acid, 75cm Long, Violet, 2/0, 37mm Needle, 1/2 Circle, RB, DCT-1
PGA283019F4P	DemeSORB, Polyglycolic Acid, 75cm Long, Violet, 3/0, 19mm Needle, 3/8 Circle, RC, DFS-2
PGA283022B0P	DemeSORB, Polyglycolic Acid, 75cm Long, Violet, 3/0, 22mm Needle, 1/2 Circle, RB, DSH-1
PGA283022B4P	DemeSORB, Polyglycolic Acid, 75cm Long, Violet, 3/0, 22mm Needle, 1/2 Circle, RC, DX-1
PGA283024F4P	DemeSORB, Polyglycolic Acid, 75cm Long, Violet, 3/0, 24mm Needle, 3/8 Circle, RC, DFS-1
PGA283026B0P	DemeSORB, Polyglycolic Acid, 75cm Long, Violet, 3/0, 26mm Needle, 1/2 Circle, RB, DSH
PGA283026B1P	DemeSORB, Polyglycolic Acid, 75cm Long, Violet, 3/0, 26mm Needle, DCT-2, 1/2 Circle, RB Heavy
PGA284019F4P	DemeSORB, Polyglycolic Acid, 75cm Long, Violet, 4/0, 19mm Needle, 3/8 Circle, RC, DFS-2
PGA284022B0P	DemeSORB, Polyglycolic Acid, 75cm Long, Violet, 4/0, 22mm Needle, 1/2 Circle, RB, DSH-1
PGA284024F4P	DemeSORB, Polyglycolic Acid, 75cm Long, Violet, 4/0, 24mm Needle, 3/8 Circle, RC, DFS-1
PGA285017B0P	DemeSORB, Polyglycolic Acid, 75cm Long, Violet, 5/0, 17mm Needle, 1/2 Circle, RB, DRB-1
PGA285019F4P	DemeSORB, Polyglycolic Acid, 75cm Long, Violet, 5/0, 19mm Needle, 3/8 Circle, RC, DFS-2
PGR1074024F4P	DemeSORB, Rapid Absorbable, 90cm Long, Colorless, 4/0, 24mm Needle, 3/8 Circle, RC, DFS-1
G1082030F4P	DemeSORB, Polyglactin 910, 90, Violet, 2/0, 30, DFSL, 3/8 Circle, RC

Polydioxanone Suture



Characteristics of DemeTech's Polydioxanone Sutures

-



Monofilament



Complete Absorption

182 - 238 Days



12 Sutures



Polydioxanone Suture

ITEM #	DESCRIPTION
PX1080036B0P	DemeDIOX, Polydioxanone, 90cm Long, Violet, 0, 36mm Needle, 1/2 Circle, RB, DCT-1
PX1082026B0P	DemeDIOX, Polydioxanone, 90cm Long, Violet, 2/0, 26mm Needle, 1/2 Circle, RB, DSH
PX1083026B0P	DemeDIOX, Polydioxanone, 90cm Long, Violet, 3/0, 26mm Needle, 1/2 Circle, RB, DSH
PX1180026B0P	DemeDIOX, Polydioxanone, 70cm Long, Violet, 0, 26mm Needle, 1/2 Circle, RB, DSH
PX1180036B4P	DemeDIOX, Polydioxanone, 70cm Long, Violet, 0, 36mm Needle, 1/2 Circle, RC, DCP-1
PX1180136B4P	DemeDIOX, Polydioxanone, 70cm Long, Violet, 1, 36mm Needle, 1/2 Circle, RC, DCP-1
PX1182026B0P	DemeDIOX, Polydioxanone, 70cm Long, Violet, 2/0, 26mm Needle, 1/2 Circle, RB, DSH
PX1182026B1P	DemeDIOX, Polydioxanone, 70cm Long, Violet, 2/0, 26mm Needle, 1/2 Circle, RB, Heavy, DCT-2
PX1182036B0P	DemeDIOX, Polydioxanone, 70cm Long, Violet, 2/0, 36mm Needle, 1/2 Circle, RB, DCT-1
PX1182036B4P	DemeDIOX, Polydioxanone, 70cm Long, Violet, 2/0, 36mm Needle, 1/2 Circle, RC, DCP-1
PX1183024F4P	DemeDIOX, Polydioxanone, 70cm Long, Violet, 3/0, 24mm Needle, 3/8 Circle, RC, DFS-1
PX1183037B0P	DemeDIOX, Polydioxanone, 70cm Long, Violet, 3/0, 37mm Needle, 1/2 Circle, RB, DCT-1
PX1184017B0P	DemeDIOX, Polydioxanone, 70cm Long, Violet, 4/0, 17mm Needle, 1/2 Circle, RB, DRB-1
PX1184022B0P	DemeDIOX, Polydioxanone, 70cm Long, Violet, 4/0, 22mm Needle, 1/2 Circle, RB, DSH-1
PX184013F13M	DemeDIOX, Polydioxanone, 45cm L, Violet, 4/0, 13mm Needle, 3/8 Circle, Precision Point RC, DP-3
PX184019F13M	DemeDIOX, Polydioxanone, 45cm Long, Violet, 4/0, 19mm Needle, 3/8 Circle Precision Point RC, DPS-2
PX185013F13M	DemeDIOX, Polydioxanone, 45cm Long, Violet, 5/0, 13mm Needle, 3/8 Circle Precision Point RC, DP-3
PX186011F0P	DemeDIOX, Polydioxanone, 45cm Long, Violet, 6/0, 11mm Needle, 3/8 Circle, RB, DBV
PX187011F0P	DemeDIOX, Polydioxanone, 45cm Long, Violet, 7/0, 11mm Needle, 3/8 Circle, RB
PX280026B0P	DemeDIOX, Polydioxanone, 75cm Long, Violet, 0, 26mm Needle, 1/2 Circle, RB, DSH
PX280026B1P	DemeDIOX, Polydioxanone, 75cm Long, Violet, 0, 26mm Needle, 1/2 Circle, RB, Heavy, DCT-2
PX280026B4P	DemeDIOX, Polydioxanone, 75cm Long, Violet, 0, 26mm Needle, 1/2 Circle, RC, DCP-2
PX280036B0P	DemeDIOX, Polydioxanone, 75cm Long, Violet, 0, 36mm Needle, 1/2 Circle, RB, DCT-1
PX280036B4P	DemeDIOX, Polydioxanone, 75cm Long, Violet, 0, 36mm Needle, 1/2 Circle, RC, DCP-1
PX280037B0P	DemeDIOX, Polydioxanone, 75cm Long, Violet, 0, 37mm Needle, 1/2 Circle, RB, DCT-1
PX280126B0P	DemeDIOX, Polydioxanone, 75cm Long, Violet, 1, 26mm Needle, 1/2 Circle, RB, DSH
PX280136B0P	DemeDIOX, Polydioxanone, 75cm Long, Violet, 1, 36mm Needle, 1/2 Circle, RB, DCT-1
PX280240B4P	DemeDIOX, Polydioxanone, 75cm Long, Violet, 2, 40mm Needle, 1/2 Circle, RC, DCP
PX282024F4P	DemeDIOX, Polydioxanone, 75cm Long, Violet, 2/0, 24mm Needle, 3/8 Circle, RC, DFS-1
PX282026B0P	DemeDIOX, Polydioxanone, 75cm Long, Violet, 2/0, 26mm Needle, 1/2 Circle, RB, DSH
PX282026B1P	DemeDIOX, Polydioxanone, 75cm Long, Violet, 2/0, 26mm Needle, 1/2 Circle, RB, Heavy, DCT-2
PX282026B4P	DemeDIOX, Polydioxanone, 75cm Long, Violet, 2/0, 26mm Needle, 1/2 Circle, RC, DCP-2
PX282026H0P	DemeDIOX, Polydioxanone, 75cm Long, Violet, 2/0, 26mm Needle, DUR-6, 5/8 Circle, RB
PX282036B0P	DemeDIOX, Polydioxanone, 75cm Long, Violet, 2/0, 36mm Needle, 1/2 Circle, RB, DCT-1
PX282036B4P	DemeDIOX, Polydioxanone, 75cm Long, Violet, 2/0, 36mm Needle, 1/2 Circle, RC, DCP-1
PX282037B0P	DemeDIOX, Polydioxanone, 75cm Long, Violet, 2/0, 37mm Needle, 1/2 Circle, RB, DCT-1
PX283017B0P	DemeDIOX, Polydioxanone, 75cm Long, Violet, 3/0, 17mm Needle, 1/2 Circle, RB, DRB-1
PX283019F4P	DemeDIOX, Polydioxanone, 75cm Long, Violet, 3/0, 19mm Needle, 3/8 Circle, RC, DFS-2
PX283022B0P	DemeDIOX, Polydioxanone, 75cm Long, Violet, 3/0, 22mm Needle, 1/2 Circle, RB, DSH-1
PX283024F4P	DemeDIOX, Polydioxanone, 75cm Long, Violet, 3/0, 24mm Needle, 3/8 Circle, RC, DFS-1
PX283026B0P	DemeDIOX, Polydioxanone, 75cm Long, Violet, 3/0, 26mm Needle, 1/2 Circle, RB, DSH
PX283026B1P	DemeDIOX, Polydioxanone, 75cm Long, Violet, 3/0, 26mm Needle, 1/2 Circle, RB, Heavy, DCT-2
PX283026B4P	DemeDIOX, Polydioxanone, 75cm Long, Violet, 3/0, 26mm Needle, 1/2 Circle, RC, DCP-2
PX283026F4P	DemeDIOX, Polydioxanone, 75cm Long, Violet, 3/0, 26mm Needle, 3/8 Circle, RC, DFS
PX283036B0P	DemeDIOX, Polydioxanone, 75cm Long, Violet, 3/0, 36mm Needle, 1/2 Circle, RB, DCT-1
PX283037B0P	DemeDIOX, Polydioxanone, 75cm Long, Violet, 3/0, 37mm Needle, 1/2 Circle, RB, DCT-1
PX284017B0P	DemeDIOX, Polydioxanone, 75cm Long, Violet, 4/0, 17mm Needle, 1/2 Circle, RB, DRB-1
PX284019F4P	DemeDIOX, Polydioxanone, 75cm Long, Violet, 4/0, 19mm Needle, 3/8 Circle, RC, DFS-2
PX284022B0P	DemeDIOX, Polydioxanone, 75cm Long, Violet, 4/0, 22mm Needle, 1/2 Circle, RB, DSH-1
PX284026B0P	DemeDIOX, Polydioxanone, 75cm Long, Violet, 4/0, 26mm Needle, 1/2 Circle, RB, DSH
PX285017B0P	DemeDIOX, Polydioxanone, 75cm Long, Violet, 5/0, 17mm Needle, 1/2 Circle, RB, DRB-1
PX285019F4P	DemeDIOX, Polydioxanone, 75cm Long, Violet, 5/0, 19mm Needle, 3/8 Circle, RC, DFS-2
PX580048B0P	DemeDIOX, Polydioxanone, 150cm Long, Violet, 0, 48mm Needle, 1/2 Circle, RB, DCTX

*Available with an Antibacterial Treatment. View Product Information on Page 26

DemeGUT™

Chronic Catgut Sutures



DemeTech's Chronic Catgut Sutures

DemeTech's Chronic Catgut suture is an absorbable, sterile, surgical suture composed of highly purified connective tissue (mostly collagen) derived from either beef or sheep intestines. DemeGUT Chronic Catgut is treated with chromium salt solution to resist body enzymes, thereby prolonging the absorption time to over 90 days. This chromic covering not only prolongs the strand's integrity, but also minimizes tissue irritation, causing less internal reaction. The percentage of collagen in the catgut suture often determines the quality of the suture. Higher percentage of collagen allows for: superior tensile strength, longer absorption times and lower reaction in vivo.

Characteristics of DemeTech's Chronic Catgut Sutures

- **DemeTech manufactures** its Chronic Catgut sutures with the highest quality collagen, maintaining a composition of 97%-98% of pure collagen.
- **DemeGUT™ Chronic Catgut sutures** are covered in chromic salt to prolong integrity, minimize tissue reactivity, and minimize the chance of adverse reactions.
- **Tensile strength** is retained for up to 14 days, and sometimes as long as 21 days.
- **Allows for smooth passage** through tissue
- **Easy knot tying** characteristics
- **DemeTech's new technology** of chromicizing strands minimizes the possibility of fraying and maximizes the potential and integrity of its Chronic Catgut.

Tensile strength
retained up to 21
days



Tensile Strength

Monofilament



Brown



Material

Complete
Absorption

56 ► 72 Days



Absorption Profile


Chromic Catgut Suture

ITEM #	DESCRIPTION
CC1124026B0P	DemeGUT, Chromic Catgut, 70cm Long, Brown, 4/0, 26mm Needle, 1/2 Circle, RB, DSH
CC124019F4P	DemeGUT, Chromic Catgut, 45cm Long, Brown, 4/0, 19mm Needle, 3/8 Circle, RC, DFS-2
CC220026B0P	DemeGUT, Chromic Catgut, 75cm Long, Brown, 0, 26mm Needle, 1/2 Circle, RB, DSH
CC222024F4P	DemeGUT, Chromic Catgut, 75cm Long, Brown, 2/0, 24mm Needle, 3/8 Circle, RC, DFS-1
CC222026B0P	DemeGUT, Chromic Catgut, 75cm Long, Brown, 2/0, 26mm Needle, 1/2 Circle, RB, DSH
CC223019F4P	DemeGUT, Chromic Catgut, 75cm Long, Brown, 3/0, 19mm Needle, 3/8 Circle, RC, DFS-2
CC224019F4P	DemeGUT, Chromic Catgut, 75cm Long, Brown, 4/0, 19mm Needle, 3/8 Circle, RC, DFS-2
CC224022B0P	DemeGUT, Chromic Catgut, 75cm Long, Brown, 4/0, 22mm Needle, 1/2 Circle, RB, DSH-1



DemeCAPRONE™

Poliglecaprone 25 Suture





DemeTech's Poliglecaprone 25 Sutures

DemeTech's fast absorbing Poliglecaprone 25 is made of a synthetic absorbable monofilament composed of glycolide and e-caprolactone copolymer (PGCL).

DemeCAPRONE is available in either violet or undyed.


Characteristics of DemeTech's Poliglecaprone 25 Sutures

- Very high tensile strength initially
- Excellent handling properties and high pliability
- Very low incidence of infection and trauma
- Smooth tissue passage
- Reliable absorption
- High in vivo strength retention



Tensile Strength

Monofilament



Material

Complete Absorption

90 ▶ 120 Days



Absorption Profile


Poliglecaprone 25 Suture

ITEM #	DESCRIPTION
MO1080030F4P	DemeCAPRONE, Poliglecaprone 25, 90cm Long, Violet, 0, 30mm Needle, 3/8 Circle, RC, DFSL
MO1080036B4P	DemeCAPRONE, Poliglecaprone 25, 90cm Long, Violet, 0, 36mm Needle, 1/2 Circle, RC, DCP-1
MO1080037B0P	DemeCAPRONE, Poliglecaprone 25, 90cm Long, Violet, 0, 37mm Needle, 1/2 Circle, RB, DCT-1
MO1080137B0P	DemeCAPRONE, Poliglecaprone 25, 90cm Long, Violet, 1, 37mm Needle, 1/2 Circle, RB
MO1080140B4P	DemeCAPRONE, Poliglecaprone 25, 90cm Long, Violet, 1, 40mm Needle, DCP, 1/2 Circle, RC
MO1080236B0P	DemeCAPRONE, Poliglecaprone 25, 90cm Long, Violet, 2, 36mm Needle, 1/2 Circle, RB, DCT-1
MO1082022F4P	DemeCAPRONE, Poliglecaprone 25, 90cm Long, Violet, 2/0, 22mm Needle, 3/8 Circle, RC
MO1082024F4P	DemeCAPRONE, Poliglecaprone 25, 90cm Long, Violet, 2/0, 24mm Needle, 3/8 Circle, RC, DFS-1
MO1082026B0P	DemeCAPRONE, Poliglecaprone 25, 90cm Long, Violet, 2/0, 26mm Needle, 1/2 Circle, RB, DSH
MO1082037B0P	DemeCAPRONE, Poliglecaprone 25, 90cm Long, Violet, 2/0, 37mm Needle, 1/2 Circle, RB, DCT-1
MO1083019F4P	DemeCAPRONE, Poliglecaprone 25, 90cm Long, Violet, 3/0, 19mm Needle, 3/8 Circle, RC, DFS-2
MO1083024F4P	DemeCAPRONE, Poliglecaprone 25, 90cm Long, Violet, 3/0, 24mm Needle, 3/8 Circle, RC, DFS-1
MO1083026B0P	DemeCAPRONE, Poliglecaprone 25, 90cm Long, Violet, 3/0, 26mm Needle, 1/2 Circle, RB
MO1083026H0P	DemeCAPRONE, Poliglecaprone 25, 90cm Long, Violet, 3/0, 26mm Needle, 5/8 Circle, RB, DUR-6
MO1083037B0P	DemeCAPRONE, Poliglecaprone 25, 90cm Long, Violet, 3/0, 37mm Needle, 1/2 Circle, RB, DCT-1
MO1084019F4P	DemeCAPRONE, Poliglecaprone 25, 90cm Long, Violet, 4/0, 19mm Needle, 3/8 Circle, RC, DFS-2
MO1182026B0P	DemeCAPRONE, Poliglecaprone 25, 70cm Long, Violet, 2/0, 26mm Needle, 1/2 Circle, RB, DSH
MO1183026B0P	DemeCAPRONE, Poliglecaprone 25, 70cm Long, Violet, 3/0, 26mm Needle, 1/2 Circle, RB, DSH
MO1184026B0P	DemeCAPRONE, Poliglecaprone 25, 70cm Long, Violet, 4/0, 26mm Needle, 1/2 Circle, RB, DSH
MO174013F13M	DemeCAPRONE Poliglecaprone 25, 45cm L, Colorless, 4/0, 13mm Needle, 3/8 Circle Precision Pt RC, DP-3
MO175013F13M	DemeCAPRONE Poliglecaprone 25, 45cm L, Colorless, 5/0, 13mm Needle, 3/8 Circle Precision Pt RC, DP-3
MO184013F13M	DemeCAPRONE Poliglecaprone 25, 45cm L, Violet, 4/0, 13mm Needle, 3/8 Circle Precision Point RC, DP-3
MO185013F13M	DemeCAPRONE Poliglecaprone 25, 45cm L, Violet, 5/0, 13mm Needle, 3/8 Circle Precision Point RC, DP-3
MO280026B1P	DemeCAPRONE, Poliglecaprone 25, 75cm Long, Violet, 0, 26mm Needle, 1/2 Circle, RB, Heavy, DCT-2
MO280030F4P	DemeCAPRONE, Poliglecaprone 25, 75cm Long, Violet, 0, 30mm Needle, 3/8 Circle, RC, DFSL
MO280236B0P	DemeCAPRONE, Poliglecaprone 25, 75cm Long, Violet, 2, 36mm Needle, 1/2 Circle, RB, DCT-1
MO282024F4P	DemeCAPRONE, Poliglecaprone 25, 75cm Long, Violet, 2/0, 24mm Needle, 3/8 Circle, RC
MO282026B0P	DemeCAPRONE, Poliglecaprone 25, 75cm Long, Violet, 2/0, 26mm Needle, 1/2 Circle, RB, DSH
MO282026B1P	DemeCAPRONE, Poliglecaprone 25, 75cm Long, Violet, 2/0, 26mm Needle, 1/2 Circle, RB, Heavy, DCT-2
MO283017B0P	DemeCAPRONE, Poliglecaprone 25, 75cm Long, Violet, 3/0, 17mm Needle, 1/2 Circle, RB, DRB-1
MO283019F4P	DemeCAPRONE, Poliglecaprone 25, 75cm Long, Violet, 3/0, 19mm Needle, 3/8 Circle, RC, DFS-2
MO283022B0P	DemeCAPRONE, Poliglecaprone 25, 75cm Long, Violet, 3/0, 22mm Needle, 1/2 Circle, RB, DSH-1
MO283024F4P	DemeCAPRONE, Poliglecaprone 25, 75cm Long, Violet, 3/0, 24mm Needle, 3/8 Circle, RC, DFS-1
MO283026B0P	DemeCAPRONE, Poliglecaprone 25, 75cm Long, Violet, 3/0, 26mm Needle, 1/2 Circle, RB, DSH
MO283026B1P	DemeCAPRONE, Poliglecaprone 25, 75cm Long, Violet, 3/0, 26mm Needle, 1/2 Circle, RB, Heavy, DCT-2
MO283026B4P	DemeCAPRONE, Poliglecaprone 25, 75cm Long, Violet, 3/0, 26mm Needle, 1/2 Circle, RC, DCP-2
MO283026F4P	DemeCAPRONE, Poliglecaprone 25, 75cm Long, Violet, 3/0, 26mm Needle, 3/8 Circle, RC, DFS
MO284017B0P	DemeCAPRONE, Poliglecaprone 25, 75cm Long, Violet, 4/0, 17mm Needle, 1/2 Circle, RB, DRB-1
MO284019F4P	DemeCAPRONE, Poliglecaprone 25, 75cm Long, Violet, 4/0, 19mm Needle, 3/8 Circle, RC, DFS-2
MO284022B0P	DemeCAPRONE, Poliglecaprone 25, 75cm Long, Violet, 4/0, 22mm Needle, 1/2 Circle, RB, DSH-1
MO284026B0P	DemeCAPRONE, Poliglecaprone 25, 75cm Long, Violet, 4/0, 26mm Needle, 1/2 Circle, RB, DSH
MO285017B0P	DemeCAPRONE, Poliglecaprone 25, 75cm Long, Violet, 5/0, 17mm Needle, 1/2 Circle, RB, DRB-1

*Available with an Antibacterial Treatment. View Product Information on Page 26


DemeCRYL™

Polyglactin 910 Suture




DemeTech's Polyglactin 910 Sutures

DemeCRYL Polyglactin 910 suture is a synthetic absorbable, sterile, surgical suture composed of copolymers made from 90% glycolide and 10% L-Lactide. The chemical formula of DemeTech's Polyglactin 910 is $(C_2H_2O_2)_m(C_3H_4O_2)_n$. Clinical trials have shown that after two weeks, approximately 75% of DemeTech's Polyglactin 910 initial strength remained. At four weeks, approximately 25% of DemeTech's Polyglactin 910 linear tensile strength remained. DemeTech's Polyglactin 910 has found to be nonantidetic and eliciting only mild tissue reactivity during the absorption process. DemeTech's Polyglactin 910 sutures are available violet and colorless.



Characteristics of Polyglactin 910 Sutures

- **Polyglactin 910 is coated** in order to create a smoother synthetic absorbable suture that will pass through tissue readily.
- **Coated Polyglactin 910 suture facilitates:**
 - Easy Tissue passage
 - Precise knot placement
 - Smooth tie down
 - A decreased tendency to irritate tissue
 - Greater tensile strength
 - A longer absorption period
- **After 3 weeks post operation**, in clinical trials, DemeTech's Polyglactin 910 sutures still maintain approximately 50% of their tensile strength.



Tensile Strength

Multifilament



Violet

Colorless



Material

Complete Absorption


56 ► 72 Days



Absorption Profile



Polyglactin Suture

ITEM #	DESCRIPTION
G1183026H0P	DemeCRYL, Polyglycolic 910, 70cm Long, Violet, 3/0, 26mm Needle, DUR-6, 5/8 Circle, RB
G273026B0P	DemeCRYL, Polyglactin 910, 75cm Long, Colorless, 3/0, 26mm Needle, 1/2 Circle, RB
G1084019F4P	DemeCRYL, Polyglactin 910, 90cm Long, Violet, 4/0, 19mm Needle, 3/8 Circle, RC, DFS-2
G1180037B0P	DemeCRYL, Polyglactin 910, 70cm Long, Violet, 0, 37mm Needle, 1/2 Circle, RB, DCT-1
G185013F13M	DemeCRYL, Polyglactin 910, 45cm Long, Violet, 5/0, 13mm Needle, 3/8 Circle Precision Point RC, DP-3
G186011F13M	DemeCRYL, Polyglactin 910, 45cm Long, Violet, 6/0, 11mm Needle, 3/8 Circle Precision Point RC, DP-1
G976011F13M	DemeCRYL, Polyglactin 910, 50cm L, Colorless, 6/0, 11mm Needle, 3/8 Circle Precision Point RC, DP-1



DemeLON™

Nylon Suture (Monofilament)



DemeTech's Nylon Sutures

DemeLON Nylon Suture is a non-absorbable, sterile surgical suture composed of the long chain aliphatic polymer Nylon 6.6. DemeLON™ Nylon sutures are available in either black or blue, dyed with FDA listed dye(s). This process is usually followed by gradual encapsulation of the suture by fibrous connective tissue.

Characteristics of DemeTech's Nylon Sutures (Monofilament)


- **Easy passage** through tissues due to its low coefficient of friction
- **Minimal** tissue reactivity
- **Excellent** elasticity
- **Exceptional** for skin closure
- **Small USP sizes** are designed specifically for micro-surgical and ophthalmic procedures

Gradual loss of
tensile strength
over time



Tensile Strength

Monofilament



Black

Blue



Material

Non
Absorbable
Suture



Absorption Profile


Nylon Suture (Monofilament)

ITEM #	DESCRIPTION
NL1063019F4P	DemeLON, Nylon Monofilament, 90cm Long, Black, 3/0, 19mm Needle, 3/8 Circle, RC, DFS-2
NL1063024F4P	DemeLON, Nylon Monofilament, 90cm Long, Black, 3/0, 24mm Needle, 3/8 Circle, RC, DFS-1
NL1064019F4P	DemeLON, Nylon Monofilament, 90cm Long, Black, 4/0, 19mm Needle, DFS-2, 3/8 Circle, RC
NL1162026B0P	DemeLON, Nylon Monofilament, 70cm Long, Black, 2/0, 26mm Needle, 1/2 Circle, RB, DSH
NL1163024F4P	DemeLON, Nylon Monofilament, 70cm Long, Black, 3/0, 24mm Needle, 3/8 Circle, RC, DFS-1
NL1164019F4P	DemeLON, Nylon Monofilament, 70cm Long, Black, 4/0, 19mm Needle, 3/8 Circle, RC, DFS-2
NL162026F4P	DemeLON, Nylon Monofilament, 45cm Long, Black, 2/0, 26mm Needle, 3/8 Circle, RC, DFS
NL163019F13M	DemeLON Nylon Monofilament, 45cm Long, Black, 3/0, 19mm Needle, 3/8 Circle Precision Point RC, DPS-2
NL163019F4P	DemeLON, Nylon Monofilament, 45cm Long, Black, 3/0, 19mm Needle, 3/8 Circle, RC, DFS-2
NL163022B4P	DemeLON, Nylon Monofilament, 45cm Long, Black, 3/0, 22mm Needle, 1/2 Circle, RC, DX-1
NL163024F4P	DemeLON, Nylon Monofilament, 45cm Long, Black, 3/0, 24mm Needle, 3/8 Circle, RC, DFS-1
NL164019F13M	DemeLON Nylon Monofilament, 45cm Long, Black, 4/0, 19mm Needle, 3/8 Circle Precision Point RC, DPS-2
NL164019F4P	DemeLON, Nylon Monofilament, 45cm Long, Black, 4/0, 19mm Needle, 3/8 Circle, RC, DFS-2
NL165013F4P	DemeLON, Nylon Monofilament, 45cm Long, Black, 5/0, 13mm Needle, 3/8 Circle, RC, DC-3
NL165019F4P	DemeLON, Nylon Monofilament, 45cm Long, Black, 5/0, 19mm Needle, 3/8 Circle, RC, DFS-2
NL192024F4P	DemeLON, Nylon Monofilament, 45cm Long, Blue, 2/0, 24mm Needle, 3/8 Circle, RC, DFS-1
NL193022B4P	DemeLON, Nylon Monofilament, 45cm Long, Blue, 3/0, 22mm Needle, 1/2 Circle, RC, DX-1
NL260148B0P	DemeLON, Nylon Monofilament, 75cm Long Black, 1, 48mm Needle, 1/2 Circle, RB
NL260240B6P	DemeLON, Nylon Monofilament, 75cm Long, Black, 2, 40mm Needle, 1/2 Circle, Taper Cutting DV-37
NL262036F4P	DemeLON, Nylon Monofilament 75cm Long, Black, 2/0, 36mm Needle, 3/8 Circle, RC, DFSL-36
NL262060A4P	DemeLON, Nylon Monofilament, 75cm Long, Black, 2/0, 60mm Needle, DKS, Straight, RC
NL262060A6P	DemeLON, Nylon Monofilament, 75cm Long, Black, 2/0, 60mm Needle, Straight Taper Cutting
NL263020B0P	DemeLON, Nylon Monofilament, 75cm Long, Black, 3/0, 20mm Needle, 1/2 Circle, RB, DSH-2
NL263024F4P	DemeLON, Nylon Monofilament, 75cm Long, Black, 3/0, 24mm Needle, 3/8 Circle, RC, DFS-1
NL263036B4P	DemeLON, Nylon Monofilament, 75cm Long, Black, 3/0, 36mm Needle, 1/2 Circle, RC, DCP-1
NL263060A4P	DemeLON, Nylon Monofilament, 75cm Long, Black, 3/0, 60mm Needle, DKS, Straight, RC
NL290030F4P	DemeLON, Nylon Monofilament, 75cm Long, Blue, 0, 30mm Needle, 3/8 Circle, RC, DFSL
NL290030F6P	DemeLON, Nylon Monofilament, 75cm Long, Blue, 0, 30mm Needle, 3/8 Circle, Taper Cutting
NL292026F4P	DemeLON, Nylon Monofilament, 75cm Long, Blue, 2/0, 26mm Needle, 3/8 Circle, RC, DFS
NL293024F4P	DemeLON, Nylon Monofilament, 75cm Long, Blue, 3/0, 24mm Needle, 3/8 Circle, RC, DFS-1
NL294019F4P	DemeLON, Nylon Monofilament, 75cm Long, Blue, 4/0, 19mm Needle, 3/8 Circle, RC, DFS-2
NL960275F4P	DemeLON, Nylon Monofilament, 50cm Long, Black, 2, 75mm Needle, 3/8 Circle, RC, DLR
NM260340B6P	DemeLON, Nylon Multifilament, 75cm Long, Black, 3, 40mm Needle, 1/2 Circle, Taper Cutting DV-37



DemeLENE™

Polypropylene Suture




DemeTech's Polypropylene Sutures

DemeLENE Polypropylene suture is a non-absorbable, sterile surgical suture composed of a synthetic linear polyolefin. DemeTech's Polypropylene is blue to enhance visibility in the body. DemeTech's Polypropylene suture elicits a minimal acute inflammatory reaction in the tissue, followed by gradual encapsulation of the suture by fibrous tissue. Unlike absorbable, and some non-absorbable sutures, DemeTech's Polypropylene is not affected or weakened by tissue enzymes. Due to its monofilament composition, Polypropylene suture resists involvement in infection. Moreover, Polypropylene has been successfully implemented in contaminated and previously infected wounds.

Characteristics of DemeTech's Polypropylene Sutures

- **Extremely high tensile strength**
- **Very flexible**
- **Due to its synthetic and monofilament characteristics,** DemeLENE™ Polypropylene elicits only minimal tissue reactivity.
- **DemeTech's Polypropylene is highly ductile.** It extends 30% prior to rupture of the thread, allowing for the post-operative swelling of tissue and avoiding possible strangulation of tissue.
- **Allows for secured knot tying**
- **Extremely smooth textile,** resulting in minimal tissue trauma



Indefinite



Tensile Strength

Monofilament



Blue

Pink



Material

Non
Absorbable
Suture



Absorption Profile

Polypropylene Suture

ITEM #	DESCRIPTION
PM1090240B0P	DemeLENE, Polypropylene, 90cm Long, Blue, 2, 40mm Needle, 1/2 Circle, RB, DCT
PM1090245B4P	DemeLENE, Polypropylene, 90cm Long, Blue, 2, 45mm Needle, 1/2 Circle, RC, DCP45
PM1094026C0P	DemeLENE, Polypropylene, 90cm Long, Blue, 4/0, 26mm Needle, 1/2 Circle, (DOUBLE) RB, DSH
PM1190037B0P	DemeLENE, Polypropylene, 70cm Long, Blue, 0, 37mm Needle, 1/2 Circle, RB, DCT-1
PM1190137B0P	DemeLENE, Polypropylene, 70cm Long, Blue, 1, 37mm Needle, 1/2 Circle, RB, DCT-1
PM192026F4P	DemeLENE, Polypropylene, 45cm Long, Blue, 2/0, 26mm Needle, 3/8 Circle, RC, DFS
PM193019F4P	DemeLENE, Polypropylene, 45cm Long, Blue, 3/0, 19mm Needle, 3/8 Circle, RC, DFS-2
PM193024F4P	DemeLENE, Polypropylene, 45cm Long, Blue, 3/0, 24mm Needle, 3/8 Circle, RC, DFS-1
PM194019F4P	DemeLENE, Polypropylene, 45cm Long, Blue, 4/0, 19mm Needle, 3/8 Circle, RC, DFS-2
PM195019F13M	DemeLENE, Polypropylene, 45cm Long, Blue, 5/0, 19mm Needle, 3/8 Circle, Precision Point RC, DPS-2
PM196013F13M	DemeLENE, Polypropylene, 45cm Long, Blue, 6/0, 13mm Needle, 3/8 Circle, Precision Point RC, DP-3
PM2110026F4P	DemeLENE, Polypropylene, 75cm Long, Pink, 0, 26mm Needle, 3/8 Circle, RC, DFS
PM2112026F4P	DemeLENE, Polypropylene, 75cm Long, Pink, 2/0, 26mm Needle, 3/8 Circle, RC, DFS
PM2113024F4P	DemeLENE, Polypropylene, 75cm Long, Pink, 3/0, 24mm Needle, 3/8 Circle, RC, DFS-1
PM2113030F4P	DemeLENE, Polypropylene, 75cm Long, Pink, 3/0, 30mm Needle, 3/8 Circle, RC, DFSL
PM2114024F4P	DemeLENE, Polypropylene, 75cm Long, Pink, 4/0, 24mm Needle, 3/8 Circle, RC, DFS-1
PM290026B0P	DemeLENE, Polypropylene, 75cm Long, Blue, 0, 26mm Needle, 1/2 Circle, RB, DSH
PM290036B0P	DemeLENE, Polypropylene, 75cm Long, Blue, 0, 36mm Needle, 1/2 Circle, RB, DCT-1
PM290037B0P	DemeLENE, Polypropylene, 75cm Long, Blue, 0, 37mm Needle, 1/2 Circle, RB, DCT-1
PM290040F4P	DemeLENE, Polypropylene, 75cm Long, Blue, 0, 40mm Needle, 3/8 Circle, RC, DFSLX
PM290040F5P	DemeLENE, Polypropylene, 75cm Long, Blue, 0, 40mm Needle, 3/8 Circle, RC, Heavy
PM290136B0P	DemeLENE, Polypropylene, 75cm Long, Blue, 1, 36mm Needle, 1/2 Circle, RB, DCT-1
PM290140B0P	DemeLENE, Polypropylene, 75cm Long, Blue, 1, 40mm Needle, 1/2 Circle, RB, DCT
PM290240B6P	DemeLENE, Polypropylene, 75cm Long, Blue, 2, 40mm Needle, 1/2 Circle, Taper Cutting DV-37
PM290340B6P	DemeLENE, Polypropylene, 75cm Long, Blue, 3, 40mm Needle, 1/2 Circle, Taper Cutting DV-37
PM292026B0P	DemeLENE, Polypropylene, 75cm Long, Blue, 2/0, 26mm Needle, 1/2 Circle, RB, DSH
PM292026B1P	DemeLENE, Polypropylene, 75cm Long, Blue, 2/0, 26mm Needle, 1/2 Circle, RB, Heavy, DCT-2
PM292026F4P	DemeLENE, Polypropylene, 75cm Long, Blue, 2/0, 26mm Needle, 3/8 Circle, RC, DFS
PM292036B0P	DemeLENE, Polypropylene, 75cm Long, Blue, 2/0, 36mm Needle, 1/2 Circle, RB, DCT-1
PM292037B0P	DemeLENE, Polypropylene, 75cm Long, Blue, 2/0, 37mm Needle, 1/2 Circle, RB, DCT-1
PM292040F4P	DemeLENE, Polypropylene, 75cm Long, Blue, 2/0, 40mm Needle, 3/8 Circle, RC, DFSLX
PM293017B0P	DemeLENE, Polypropylene, 75cm Long, Blue, 3/0, 17mm Needle, 1/2 Circle, RB, DRB-1
PM293019F4P	DemeLENE, Polypropylene, 75cm Long, Blue, 3/0, 19mm Needle, 3/8 Circle, RC, DFS-2
PM293024F4P	DemeLENE, Polypropylene, 75cm Long, Blue, 3/0, 24mm Needle, 3/8 Circle, RC, DFS-1
PM293026B0P	DemeLENE, Polypropylene, 75cm Long, Blue, 3/0, 26mm Needle, 1/2 Circle, RB, DSH
PM293026B1P	DemeLENE, Polypropylene, 75cm Long, Blue, 3/0, 26mm Needle, 1/2 Circle, RB, Heavy, DCT-2
PM294019B0P	DemeLENE, Polypropylene, 75cm Long, Blue, 4/0, 19mm Needle, 1/2 Circle, RB
PM3X990275F4P	DemeLENE, Polypropylene, 3x50cm Long, Blue, 2, 75mm Needle, 3/8 Circle, RC, DLR
PM590265B0P	DemeLENE, Polypropylene, 150cm Long, Blue, 2, 65mm Needle, 1/2 Circle, RB, DTP-1


“No technical problems here with the sutures!”

Pam B., Ontario - CANADA

Using PDO, PGCL, PGA, Polypro

DemeSILK™

Silk Suture



DemeTech's Silk Sutures

DemeSILK Silk suture is a non-absorbable, sterile, surgical suture composed of an organic protein called fibroin. DemeTech's silk sutures are processed to remove the natural waxes and gums from its origin. As Silk is a natural fiber, it does elicit an acute inflammatory reaction, which is followed by gradual encapsulation of the suture by fibrous connective tissue (occurring usually in 14-21 days). While DemeTech's silk sutures are not considered absorbable sutures, progressive degradation of its fibers results in the gradual loss of tensile strength within the first year. DemeTech's silk sutures are indicated for use in general soft tissues approximation and/or ligation.

Characteristics of DemeTech's Silk Sutures

- Raw silk is produced by a silkworm larva's cocoon. The natural waxes and gums are removed prior to the braided process, allowing for a tighter, more compact braid.
- DemeTech's Silk suture is dyed black for easier visibility in tissue.
- DemeSILK Silk suture is non-absorbable and multifilament in nature.
- Although DemeTech classifies its silk suture as a non-absorbable suture, it loses most of its tensile strength in one year, and is completely absorbed in 2 years. Thus, DemeTech's Silk behaves like a very slow absorbing suture.
- DemeTech coats its silk in silicon to limit absorption, increase ductility, and minimize tissue reactivity.

Gradual loss of
tensile strength
over time



Tensile Strength

Multifilament



Black



Material

Non
Absorbable
Suture



Absorption Profile

Silk Suture

ITEM #	DESCRIPTION
SK162026F4P	DemeSILK, Silk, 45cm Long, Black, 2/0, 26mm Needle, 3/8 Circle, RC, DFS
SK163024F4P	DemeSILK, Silk, 45cm Long, Black, 3/0, 24mm Needle, 3/8 Circle, RC, DFS-1
SK164019F4P	DemeSILK, Silk, 45cm Long, Black, 4/0, 19mm Needle, 3/8 Circle, RC, DFS-2
SK260025B4P	DemeSILK, Silk, 75cm Long, Black, 0, 25mm Needle, 1/2 Circle, RC
SK262017F0P	DemeSILK, Silk, 75cm Long, Black, 2/0, 17mm Needle, 3/8 Circle, RB, DBB
SK2X5600W	DemeSILK, Silk, 2x150cm Long, Black, 0mm Needle, No Needle, No Needle



Antibacterial Poliglecaprone 25 Suture

ITEM #	DESCRIPTION
MO1080030F4AB	DemeCAPRONE, Poliglecaprone 25 Antibacterial, 90cm Long, Violet, 0, 30mm Needle, 3/8 Circle RC, DFSL
MO1080036B4AB	DemeCAPRONE Poliglecaprone 25 Antibacterial, 90cm Long, Violet, 0, 36mm Needle, 1/2 Circle RC, DCP-1
MO1080037B0AB	DemeCAPRONE, Poliglecaprone 25 Antibacterial, 90cm L, Violet, 0, 37mm Needle, 1/2 Circle, RB, DCT-1
MO1082024F4AB	DemeCAPRONE, Poliglecaprone 25 Antibacterial, 90cm L, Violet, 2/0, 24mm Needle, 3/8 Circle RC, DFS-1
MO1082037B0AB	DemeCAPRONE, Poliglecaprone 25 Antibacterial, 90cm L, Violet, 2/0, 37mm Needle, 1/2 Circle RB, DCT-1
MO1083019F4AB	DemeCAPRONE, Poliglecaprone 25 Antibacterial, 90cm L, Violet, 3/0, 19mm Needle, 3/8 Circle RC, DFS-2
MO1083024F4AB	DemeCAPRONE, Poliglecaprone 25 Antibacterial, 90cm L, Violet, 3/0, 24mm Needle, 3/8 Circle RC, DFS-1
MO1084019F4AB	DemeCAPRONE, Poliglecaprone 25 Antibacterial, 90cm L, Violet, 4/0, 19mm Needle, 3/8 Circle RC, DFS-2
MO282026B0AB	DemeCAPRONE, Poliglecaprone 25 Antibacterial, 75cm L, Violet, 2/0, 26mm Needle, 1/2 Circle, RB, DSH
MO283024F4AB	DemeCAPRONE, Poliglecaprone 25 Antibacterial, 75cm L, Violet, 3/0, 24mm Needle, 3/8 Circle RC, DFS-1
MO283026B0AB	DemeCAPRONE, Poliglecaprone 25 Antibacterial, 75cm L, Violet, 3/0, 26mm Needle, 1/2 Circle, RB, DSH
MO284019F4AB	DemeCAPRONE, Poliglecaprone 25 Antibacterial, 75cm L, Violet, 4/0, 19mm Needle, 3/8 Circle RC, DFS-2
MO284022B0AB	DemeCAPRONE, Poliglecaprone 25 Antibacterial, 75cm L, Violet, 4/0, 22mm Needle, 1/2 Circle RB, DSH-1
MO284026B0AB	DemeCAPRONE, Poliglecaprone 25 Antibacterial, 75cm L, Violet, 4/0, 26mm Needle, 1/2 Circle, RB, DSH

Antibacterial Polydioxanone Suture

ITEM #	DESCRIPTION
PX1080036B0AB	DemeDIOX, Polydioxanone Antibacterial, 90cm Long, Violet, 0, 36mm Needle, 1/2 Circle, RB, DCT-1
PX280026B4AB	DemeDIOX, Polydioxanone Antibacterial, 75cm Long, Violet, 0, 26mm Needle, 1/2 Circle, RC, DCP-1
PX280036B0AB	DemeDIOX, Polydioxanone Antibacterial, 75cm Long, Violet, 0, 36mm Needle, 1/2 Circle, RB, DCT-1
PX280036B4AB	DemeDIOX, Polydioxanone Antibacterial, 75cm Long, Violet, 0, 36mm Needle, 1/2 Circle, RC, DCP-1
PX282024F4AB	DemeDIOX, Polydioxanone Antibacterial, 75cm Long, Violet, 2/0, 24mm Needle, 3/8 Circle, RC, DFS-1
PX282026B0AB	DemeDIOX, Polydioxanone Antibacterial, 75cm Long, Violet, 2/0, 26mm Needle, 1/2 Circle, RB, DSH
PX282026B4AB	DemeDIOX, Polydioxanone Antibacterial, 75cm Long, Violet, 2/0, 26mm Needle, 1/2 Circle, RC, DCP-2
PX282036B0AB	DemeDIOX, Polydioxanone Antibacterial, 75cm Long, Violet, 2/0, 36mm Needle, 1/2 Circle, RB, DCT-1
PX282036B4AB	DemeDIOX, Polydioxanone Antibacterial, 75cm Long, Violet, 2/0, 36mm Needle, 1/2 Circle, RC, DCP-1
PX283019F4AB	DemeDIOX, Polydioxanone Antibacterial, 75cm Long, Violet, 3/0, 19mm Needle, 3/8 Circle, RC, DFS-2
PX283022B0AB	DemeDIOX, Polydioxanone Antibacterial, 75cm Long, Violet, 3/0, 22mm Needle, 1/2 Circle, RB, DSH-1
PX283024F4AB	DemeDIOX, Polydioxanone Antibacterial, 75cm Long, Violet, 3/0, 24mm Needle, 3/8 Circle, RC, DFS-1
PX283026B0AB	DemeDIOX, Polydioxanone Antibacterial, 75cm Long, Violet, 3/0, 26mm Needle, 1/2 Circle, RB, DSH
PX283026B4AB	DemeDIOX, Polydioxanone Antibacterial, 75cm Long, Violet, 3/0, 26mm Needle, 1/2 Circle, RC, DCP-2
PX284019F4AB	DemeDIOX, Polydioxanone Antibacterial, 75cm Long, Violet, 4/0, 19mm Needle, 3/8 Circle, RC, DFS-2
PX284022B0AB	DemeDIOX, Polydioxanone Antibacterial, 75cm Long, Violet, 4/0, 22mm Needle, 1/2 Circle, RB, DSH-1



DemeWAX™

Sterile BoneWax 2.5gm



DemeWAX is a sterile wax used to help mechanically control bleeding from Bone surfaces during surgical procedures.


- Can be used as a local hemostat on bone surfaces.
- Non Absorbable and does not elicit biochemical reactions.
- DemeWax is composed of a sterile mixture of beeswax paraffin, and isopropyl palmitate.

DemeTech's BoneWax

ITEM #	DESCRIPTION
BW2.5	Sterile Bonewax 2.5g

DemeLENE MESH™

Polypropylene Surgical Mesh



DemeTech's Surgical Mesh

ITEM #	DESCRIPTION
DPM15X15	DemeLENE Polypropylene Surgical Mesh, Nonabsorbable, 15cm x 15cm (6" x 6")
DPM30X30	DemeLENE Polypropylene Surgical Mesh, Nonabsorbable, 30cm x 30cm (12" x 12")
DPM6X11	DemeLENE Polypropylene Surgical Mesh, Nonabsorbable, 6cm x 11cm (2.3" x 4.3")



Movora

Local contact: www.movora.com/contact-us

© 2025 Movora

Subject to technical changes.

DemeSORB, DemeDIOX, DemeGUT, DemeCAPRONE,
DemeLON, DemeLENE, and DemeWAX are trademarks
of DemeTECH, used with permission.

All rights reserved.

MMKT-PC-0032-A