

Crescentic Osteotomy Blades

X-RAY TEMPLATE

100% Magnification



8 mm

0.14mm per degree
of correction




10 mm

0.17mm per degree
of correction




12 mm

0.21mm per degree
of correction




15 mm

0.28mm per degree
of correction



18 mm

0.32mm per degree
of correction



21 mm

0.35mm per degree
of correction



24 mm

0.42mm per degree
of correction



27 mm

0.49mm per degree
of correction



30 mm

0.53mm per degree
of correction



33 mm

0.58mm per degree
of correction

NEW GENERATION DEVICES

65 Harristown Rd., Glen Rock, NJ 07452 800-797-8587 201-891-5615

www.ngdvet.com