

Versiv TPLO

Our TPLO plates were designed leveraging the R&D excellence within the Movora brands with features promoting optimal vascularity and accelerated healing.

FEATURES:

- Curved distal portion allows for placement freedom
- Hohmann retractor slot for minimally invasive incisions
- Proximal screw holes are parallel to the osteotomy providing better stability and resistance to pull out. The most proximal hole is angled down 3°
- Minimal contact v-shaped undersurface for vascularity and accelerated healing
- Divot for reduction forceps
- Use with Universal Drill Guide



Standard 3.5 mm shown

BROAD SIZE RANGE FROM TOY TO JUMBO PLUS

Available sizes:

1.5 Toy | 2.0 | 2.4 | 2.7 | 2.7 Broad | 3.5 Small | 3.5 | 3.5 Broad | 3.5 Jumbo | 4.0/5.0 Jumbo Plus

Use with Universal Drill Guide.
Accepts the following screws:

- ⊗⊗ Double Threaded Locking (SLH ST or HLH ST)
- ⊗⊗ Cortical (S ST or ST)
- ⊗ TorqAlign (TA SLH ST or TA S ST)

Versiv TPLO

Screw Hole Sequence

- 1 CORTICAL SCREW**
Partially Tightened
- 2 LOCKING SCREW**
Fully Tightened
- 3 LOCKING SCREW**
Fully Tightened
- 4 CORTICAL SCREW**
*Compression Screw
Fully Tightened*
- 5 LOCKING OR
CORTICAL SCREW**
*For Additional Compression
Fully Tightened*
- 6 LOCKING SCREW**
Fully Tightened
- 7 LOCKING SCREW**
Fully Tightened

